

The World's first portable direct LDL cholesterol test

# CardioChek PA portable whole blood analyser



**1 device. 6 different parameters. 10 different combinations.**

## Lipid Panel:

- Total Cholesterol
- Triglycerides
- HDL Cholesterol
- LDL Cholesterol (calculated)

## Metabolic Chemistry Panel:

- HDL Cholesterol
- Glucose
- Triglycerides

## CVD Panel:

- Total Cholesterol
- Glucose
- HDL Cholesterol

## General Screening Panel:

- Total Cholesterol
- HDL Cholesterol

All above are available as single test panels.

- **Innovative**

Runs on portable hand-held device

- **Fast**

Results in one to two minutes

- **Flexible**

Suitable for multiple clinical applications

- **Accurate**

CE marked/FDA cleared/CLIA waived

- **Convenient**

Test at point-of-care, no fasting, no refrigeration

- **Cost-Effective**

No expensive lab analysis or follow-up consultation delay

- **Fully automated**

Simply pipette a single sample onto strip

- **Memory log**

Stores up to 200 results – date and time stamped.

**And all in the time it takes you to do a BP test!**



Call: +44 (0) 8700 349 606 Visit: [www.bhr.co.uk](http://www.bhr.co.uk) Write: [info@bhr.co.uk](mailto:info@bhr.co.uk)

BHR Pharmaceuticals Limited, Centenary Business Centre, Hammond Close, Nuneaton, Warwickshire CV11 6RY

# CardioChek PA analyser

## technical specification

|                         |   |                              |  |
|-------------------------|---|------------------------------|--|
| <b>measuring theory</b> | reflectance photometry  | <b>measuring temperature</b> | 18-35°C  |
| <b>test range</b>       | cholesterol: 2.59-10.36mmol/L<br>hdl: 0.39-2.59mmol/L<br>ldl: 1.29-5.18mmol/L<br>triglycerides: 0.57-5.65mmol/L<br>glucose: 1.11-3.33mmol/L | <b>measuring method</b>      | spectrophotometry                                |
|                         |   | <b>measuring time</b>        | 1.6 minutes (full lipid)<br>45 seconds (glucose) |
|                         |   | <b>sample</b>                | whole blood. 15-40µl                             |

## performing a test

take a 15-40µl fingerprick sample ... dispense sample onto strip ..... read result ... simple!



## available tests

| code | test   | contents       |
|------|--|----------------|
|      | BIOS002 <b>glucose</b> (requires 15µl pipette)   | 25 test panels |
|      | BIOS003 <b>total cholesterol</b> (requires 15µl pipette)   | 25 test panels |
|      | BIOS004 <b>hdl cholesterol</b> (requires 15µl pipette)   | 25 test panels |
|      | BIOS005 <b>triglyceride</b> (requires 15µl pipette)  | 25 test panels |
|      | BIOS006 <b>ketone</b> (requires 15µl pipette)  | 25 test panels |
|      | BIOS016 <b>direct ldl</b> (requires 15µl pipette)  | 25 test panels |
|      | BIOS002D <b>total cholesterol &amp; glucose duo</b> (requires 30µl pipette)                              | 25 test panels |
|      | BIOS010 <b>total cholesterol &amp; hdl duo</b> (requires 30µl pipette)                                   | 25 test panels |
|      | BIOS019 <b>total cholesterol , hdl &amp; glucose triple</b> (requires 40µl pipette)                      | 15 test panels |
|      | BIOS020 <b>hdl, triglyceride &amp; glucose (metabolic chemistry) triple</b> (requires 40µl pipette)      | 15 test panels |
|      | BIOP102 <b>lipid panel (total cholesterol, hdl, ldl &amp; triglyceride) quad</b> (requires 40µl pipette) | 15 test panels |

## monitors

| code    | item                                | contents |
|---------|-------------------------------------|----------|
| BIOP101 | <b>cardiochek pa blood analyser</b> | 1 kit    |
| BIOP104 | <b>cardiochek blood analyser</b>    | 1 kit    |



## related items

| code    | item  | contents       |
|---------|---|----------------|
| BIOS011 | <b>roche softclix pro lancing pen</b>         | 1 device       |
| BIOS012 | <b>roche softclix pro lancets</b>             | 200 lancets    |
| BIOS009 | <b>safetec pipette 15µl</b>                   | 25 pipettes    |
| BIOS013 | <b>safetec pipette 30µl</b>                   | 25 pipettes    |
| BIOP103 | <b>safetec pipette 40µl</b>                   | 16 pipettes    |
| BIOS015 | <b>drummond tube &amp; plunger 15µl</b>       | 25 sets        |
| BIOS018 | <b>drummond tube &amp; plunger 30µl</b>       | 25 sets        |
| BIOP108 | <b>drummond tube &amp; plunger 40µl</b>       | 25 sets        |
| BIOS007 | <b>pts multichemistry control solution</b>    | 1l & 1h vials  |
| BIOS017 | <b>pts hdl-cholesterol control solution</b>   | 1l & 1h vials  |
| BIOP014 | <b>medi-wipes</b>                             | 100 wipes      |
| BIOP105 | <b>cardiochek printer</b>                     | 1 kit          |
| BIOP106 | <b>cardiochek printer communication cable</b> | 1 cable        |
| BIOP107 | <b>cardiochek printer label roll</b>          | 5 x 500 labels |

